***FACE SHIELD - V01***

## Technical Data Sheet

### Face Shield - V01

**Composition:** Polycarbonate Face Shield 1mm, PUR Foam 35mm and Headband 3mm.

---

**Description:** Polycarbonate Face Shield with Protective Foam and Adjustable Headband.  
Can be used with or without glasses and mask.

---

**Materials:**

- Crystal Polycarbonate 1mm
- Grey PUR Foam 35mm
- Black Headband 3mm

---

**Characteristics:**

- Perfect Finishing
- Comfortable
- Safe
- Reusable
- Durable

---

**Applications:**

- Hospitals
- Pharmacy
- Customer Service
- Firefighters
- Manufacturers
- Nursing Home
- Transporters

## Technical Data Sheet

### Face Shield - V01

#### Placement and adjustment:

- Wash your hands with soap and water drying them afterwards
- Place the Face Shield over your face and eyes without touching the inside, as well as ensuring that the foam sits comfortably on your forehead, slightly above your eyes
- Adjust the headband if necessary

---

#### Removal and disinfection:

- Wash your hands with soap and water, drying them afterwards
- Remove the Face Shield, manipulating only the headband. The front, inside and sides must never be touched.
- Place the Face Shield in a recipient and wash with water and a neutral soap, disinfecting afterwards with 70% alcohol. Leave to air dry
- Wash your hands with soap and water, drying them afterwards

---

#### Important Note:

Ideally gloves should be used when removing, cleaning and disinfecting.

In case you don't use gloves and your hands have been contaminated during removal of the Face Shield, immediately wash your hands or use a alcohol based disinfectant for your hands.

Always wash your hands at the end of the complete process.

Before use, please remove the film on the both sides of the Face Shield

## Technical Data Sheet

### Face Shield - V01

#### Polycarbonate Specifications:

- Density: 1,2 gr/cm<sup>3</sup> ASTM D792
- Tensile Modulus of Elasticity: 2300 MPa ASTM D638
- Service Temperature Range: -40°C a 120°C.
- Hardness: 125R
- Light Transmission: 89%
- Flexural Strength: 100 MPa
- Coefficient of Linear Thermal Expansion:  $6.5 \times 10^{-5}$  cm/cm °C
- Tensile Strength at Yield: 65 MPa ASTM D638
- Elongation at Break: >90% ASTM D638
- Impact Falling Weight: 158J ISO6603/1
- Thermal Conductivity: 0.21 W/mK ASTM C177

#### Headband Specifications:

- Length: 44cm
- Thickness: 3mm
- Composition: 34,5% polyester e 65,5% rubber

#### Foam Specifications:

- Density: 33Kg/m<sup>3</sup> DIN EN ISO 845
- Compressive Strength, 40% 2,5-4,0 kPa DIN EN ISO 3386-1
- Resistance to Traction 180kPa DIN EN ISO 1798
- Elongation at Break: 220% DIN EN ISO 1798
- Porosity: 42-50 PPI
- Temperature Range: -40° a 120°C

#### Storage conditions:

Face Protection Visors should be stored in a cool, dry place and not be exposed to extreme heat.

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